

Deluxe Folding Walker, Two Button with 5" Wheels & Tennis Ball Glides

10210TB-2 Adult, 2/cs

FEATURES:

- Comes with 5" wheels and Tennis Ball glides. (item #10119 with additional glides, 10107 and 10107C).
- Comes with wheels. For replacement wheels, refer to item # 10109.
- Easy push-button mechanisms may be operated by fingers, palm or side of hand.
- Each side operates independently to allow easy movement through narrow spaces and greater stability while standing.
- Unique U-shape frame design creates greater clearance.
- The patient is able to bring the walker closer to assist in standing.
- Comes with rear tennis ball glides.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Vinyl contoured handgrip for added comfort.
- Limited Lifetime Warranty.



10210TB-2



TENNIS GLIDES ENABLE WALKER TO SMOOTHLY GLIDE ALONG FLOOR SURFACES!

New & Improved Bariatric (Wider & Deeper Frame Design) Aluminum Folding Walker, Two Button

BARIATRIC
500 LB
WEIGHT CAPACITY



THIS ITEM IS AVAILABLE IN A 1-PACK!

- 10220-1 Adult, 1/cs
- 10220-2 Adult, 2/cs
- 10223-2 Junior, 2/cs
- 10220-2WW Adult, with Heavy Duty Wheels, 2/cs

FEATURES:

- For Bariatric use.
- Steel legs and side braces.
- Wider and deeper frame design accommodates individuals up to 500 lbs.
- Vinyl contoured hand grip.
- Item # 10220-2WW comes with wheels. For replacement wheels, refer to item # 10118SV.
- Precision design provides additional strength while adding minimal weight.
- Limited Lifetime Warranty.



10220-2WW



10220-1

TWO BUTTON BARIATRIC WALKER WITH WHEELS FEATURE MATCHING SILVER VEIN "U" BRACE AND SIDE FRAMES & LEGS!

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10210TB-2	17"	19"	25.5"	18"	39"	32"	6.9 lbs.	15.8 lbs.	350 lbs.
10220-1	20"	23"	24"	21"	39"	32"	7.5 lbs.	8.8 lbs.	500 lbs.
10220-2	20"	23"	24"	21"	39"	32"	7.5 lbs.	17 lbs.	500 lbs.
10223-2	20"	23"	24"	20"	32"	25"	7.3 lbs.	16.6 lbs.	500 lbs.
10220-2WW	20"	23"	27"	24"	39"	32"	11 lbs.	24.2 lbs.	500 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.