



300 LB
WEIGHT CAPACITY

Med Aire Alternating Pressure Pump and Pad System

14001E	78" (L) x 34" (W) x 2.5" (H), Variable Pressure Pump and Pad, (Has replaced 14001D), 1/cs
14001E-220	78" (L) x 34" (W) x 2.5" (H), Variable Pressure Pump (220V) and Pad, (Has replaced 14001D), 1/cs
14001EFD	78" (L) x 34" (W) x 2.5" (H), Variable Pressure Pump and Deluxe Pad with End Flaps, 1/cs
14002E	78" (L) x 34" (W) x 2.5" (H), Fixed Pressure Pump and Deluxe Pad, (Has replaced 14002D), 1/cs
14002E-220	78" (L) x 34" (W) x 2.5" (H), Fixed Pressure Pump (220V) and Deluxe Pad, (Has replaced 14002D), 1/cs
14003	78" (L) x 34" (W) x 2.5" (H), Standard Pad Only, 1/ea
14003-EF	78" (L) x 34" (W) x 2.5" (H), Deluxe Pad with End Flaps, 1/ea

HCPCS code: E0181

FEATURES:

- 130 individual bubble cells for maximum comfort and effectiveness.
- Quiet pump technology alternately inflates and deflates the air cells.
- **4 LPM (Litres per Minute) pump produces consistent air flow and pressure.**
- 5-Minute cycle time.
- Easy to change filter.
- Variable Pressure setting on Deluxe Pump allows comfort setting for maximum results.
- Bubble pad design provides superior therapeutic treatment.
- Built in bracket for mounting easily to bed.
- Pad is made of long lasting durable heavy gauge premium vinyl.
- 1 Year Warranty on the pump.
- 90 day Warranty on mattress.



14005E



14006E

REPLACEMENT PARTS & ACCESSORIES

14001T	Beige Tubing for 14001 and 14001E, 1/ea
14002T	Clear Tubing for 14002 and 14002E, 1/ea
14005E	Variable Pressure Pump Only for 14001E, (Has replaced 14005D) 1/ea
14006E	Fixed Pressure Pump Only for 14002E, (Has replaced 14006D), 1/ea
14005EF	Filter for 14005E, 14006E, 1/ea

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.