

# lifetime limited warranty

Your Drive brand product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser.

This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

In the event of a defect covered by this warranty, we will, at our option, repair or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



*Wenzel<sup>ite</sup>Re/hab*

© 2006 Medical Depot, Inc. All rights reserved.

Drive is a trademark of Medical Depot, Inc.

Port Washington N.Y. 11050 USA Made in China

dolphin bathing chair base

Item # DO 2000

*Wenzel<sup>ite</sup>Re/hab*

[www.drivemedical.com](http://www.drivemedical.com)

drive  
MEDICAL DESIGN & MANUFACTURING

---

## recommended use

---

Height adjustable base easily attaches to the Dolphin Chair. It allows for easy transfer and can be used in either the bathtub or shower.

---

## assembly

---

1. Remove rubber plugs from Dolphin chair frame (A).
2. Attach base to chair using four aluminum screws supplied (B).

---

## adjustments

---

1. Swing open to use or fold to close.
2. For height adjustment, raise or lower blue U-frame to desired height. Secure position by insertnig collar-pins (C)

---

## caution

---

Adult supervision is required at all times. The straps and supports must be carefully adjusted for the comfort and security of the child. Child must be strapped in at all times.



Item # DO 2010, Dolphin Tub Stand



Item # DO 2000,  
Dolphin Bathing Chair



Item # DO 2000 & DO 2010  
Dolphin Bath Chair  
With Tub Stand