



- ✓ Textured, non-slip grip bag
- ✓ Secure-fitting mask
- ✓ Flexible, 360 degree swivel



## Manual Resuscitators

Convenience, reliability and value





# Introducing Medline Manual Resuscitators

All the best features at a price that fits your budget

- Bag is textured for an enhanced, non-slip grip.
- Premium, soft cushion mask allows for an effective seal on the patient's face.
- Lightweight, translucent bag body and components allow for visualization of vomitus or occlusions in the system.
- The patient port swivels 360 degrees, allowing the clinician freedom to move around the patient without placing excessive torque or stress on the endotracheal tube or mask.
- The body of the bag swivels 360 degrees, allowing for multiple clinician hand positions, helping to reduce clinician fatigue and providing greater patient comfort.
- Bright red, universal O<sub>2</sub> tubing connector provides easy confirmation that the resuscitation bag is attached to the flowmeter.
- Infant, pediatric and adult sizes available.

## Resuscitator Configurations

Item #	Bag Size			Mask	Reservoir Style		PEEP Valve	Case Qty.
	Adult	Pediatric	Infant		Bag	Tubing		
CPRM1116	•			•	•			6
CPRM1126	•			•		•		6
CPRM1116P	•			•	•		•	6
CPRM2216		•		•	•			6
CPRM2226		•		•		•		6
CPRM3312			•	•	•			12
CPRM3322			•	•		•		12

## Fittings and Accessories



- **CPRMP19** – 19mm Adjustable PEEP Valve
- **CPRMP30** – 30mm Adjustable PEEP Valve with 19mm Adapter
- **CPRMCO2** – CO<sub>2</sub> Indicator



**Medline Industries, Inc.**  
One Medline Place  
Mundelein, IL 60060

**Medline United States**  
**1-800-MEDLINE (633-5463)**  
www.medline.com | info@medline.com

**Medline Canada**  
**1-800-396-6996**  
www.medline.ca | canada@medline.com

**Medline México**  
**01-800-831-0898**  
www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. We reserve the right to correct any errors that may occur within this brochure.  
©2011 Medline Industries Inc. Medline is a registered trademark of Medline Industries, Inc.