

DURABLE MEDICAL EQUIPMENT



Effective Programs, Quality Products,
and Superior Customer Service

Safe Bathrooms

Bath Benches

Bath Accessories

Transfer Benches

Elevated Toilet Seats

Grab Bars

Safe-T-Pole

Commodes

Making Your Bathroom Safer

To help provide security and stability, use the guide below to understand where bath safety products can be located in and around the bathroom.



A. Handheld Shower

G98035

Hand held showers should be used for anyone who bathes with a shower chair or transfer bench. The 86" (218 cm) long hose allows a person to remain seated while bathing.

B. Elevated Toilet Seat

G30270A

When bending is difficult, using the toilet can present a potential fall risk. Using a raised toilet seat will reduce the degree of required bending, thereby reducing fall risk.

C. Bathroom Grab Bars and Assist Bars

Install grab bars in your bathroom to offer support in areas where falls are most likely to occur. These areas include anywhere that might have a wet surface or anywhere people tend to bend frequently. A grab bar's primary use is to steady the user and may not support the weight of a falling individual.

G98006

Bathtub Assist Bar – Clamps directly onto the tub wall. The assist bar helps make entry and exit safer.

Wall-Mounted Grab Bars (not shown) – Mounts directly to your wall to provide an extra grip when standing or sitting in the tub.

D. Transfer Bench

MDS86952

When getting in and out of the tub becomes a difficult and potentially dangerous task, transfer benches are a careful alternative. The person bathing can sit down on the outside of the bathtub, swing their legs over the wall, and then easily slide across the seat while holding onto a hand rail mounted on the opposite side.

Bathroom Safety and Mobility Questions—Answered by Medline



Health conditions can create safety challenges in the bathroom and impair your mobility. These answers can point you towards products that improve bathroom safety and enhance your mobility.

How can I make it easier to get into and out of my bathtub?

A **transfer bench** straddles the side of the bathtub, allowing you to sit down outside the tub and then move over the surface of the bench into the bathtub. You may also want to consider installing grab bars on the walls of the shower/tub area.

How can I prevent falling when I am bathing or showering?

A **bath bench** can be placed inside your bathtub or shower. You can then safely sit on the bench while you wash your body. You may also want to consider installing grab bars on the walls and using a hand-held shower attachment to make showering easier.

How can I make it easier to get up from a toilet seat?

A **raised toilet seat** can easily be installed on top of your toilet seat. The advantage of a raised toilet seat is that it provides additional height, decreasing the distance you have to sit down or stand up.

How can I keep myself steady and prevent falling when I am standing in the shower?

Consider installing **grab bars** on the walls. Grab bars provide safety and support to help you maintain balance.

How can I keep the cast on my arm dry while I am bathing or showering?

A **cast protector** is a plastic sleeve that can be pulled over your hand and forearm, shielding your cast from the water.

How can I take care of toileting closer to my bed?

A **bedside commode** can be placed next to your bed, making it easy and accessible anytime. You can also use commode liners inside the commode bucket to make cleanup easier and sanitary.

What kind of walking aid would allow me to safely walk along outdoors, but sit down if I become tired?

A **rollator** has four wheels, handbrakes, and a padded seat on which you can sit. A rollator allows you to safely walk along, but sit on the padded seat if you become tired.

Should I use a standard wheelchair (large wheels in rear) or a transport chair (small wheels in rear)?

A **standard wheelchair** allows you to self-propel yourself without the aid of a caregiver. Although a **transport chair** requires a caregiver to push the chair, a transport chair is smaller and lighter than a standard wheelchair—making a transport chair easier to transport in an automobile. Some people choose to own both a standard wheelchair and a transport chair, so they have the best of both worlds.

What kind of walking aid would allow me to be pushed (transported) part of the time and walk along part of the time?

A **combination transport chair/rollator** is a convenient device which converts in seconds between the modes of transport chair and rollator. You can walk along and get exercise, sit and rest when you are tired, or be pushed along with the device in transport chair mode.

Bath Benches

- Bath benches help patients who have difficulty sitting in a standard bath tub or standing in the shower
- Benches with a back give extra support and comfort
- Stools without a back provide easy access and maneuverability when bathing
- Non-slip, rubber feet resist slipping on the bathtub or shower floor
- Legs are height-adjustable for a proper fit



G98027

Premium Padded Bath Bench

Extra thick padding on seat and backrest provides additional comfort for long term use. Bench features an extra deep seat 16.5" (42 cm), for greater stability and comfort when bathing.



MDS89755R

Standard Composite Plastic Bath Bench

Durable and rustproof all-plastic construction. Comes with an extra deep and wide seat for user comfort. Removable arms help aid patient balance, but can be removed for easier transfers.



G30402

Basic Bath Bench

Textured, blow-molded polyethylene seats have built-in hand holds to provide extra safety. Angled legs distribute weight over a large area to increase patient stability. Durable anodized aluminum frame is rust resistant.



MDS89745AXW

Premium Bariatric Bath Bench 550 lb Capacity

Dual cross brace with durable welds provide extra stability and support. Bench has a blow molded seat with high-strength, aluminum frame to prevent rusting.



MDS89745RA*

Standard Bath Bench with Arms

Ideal for patients that need extra support when getting on and off the bench. Padded hand-grip for comfort also reduces slipping.



MDS89745R

Basic Bath Bench

High Density, polyethylene seat and back with built-in hand grips. Angled legs distribute weight over a large area to increase patient security. Durable, anodized aluminum frame is rust resistant.

GUARDIAN® Bath Bench Accessories

- Additional items to go with your bath bench to make bathing easier



CUR200ALL



G98035



G01212

Cast Protectors

Ideal for keeping casts, bandages or wounds dry while showering. Specially designed opening stretches to fit over the cast and then returns to its original shape to fit snugly against the arm.

Item Number	Description	Pkg
CUR200AAA	Adult Arm	2/box, 6 bx/cs
CUR200ALL	Adult Leg	2/box, 6 bx/cs
CUR100KAA	Child Arm	2/box, 6 bx/cs
CUR100KLL	Child Leg	2/box, 6 bx/cs

Ergonomic Hand Held Shower

Extra long 86" (218 cm) hose allows user to bathe comfortably while seated in shower. Hand held shower easily adapts to existing shower arm, delivering standard or pulsating spray. Diverter valve allows you to use the existing shower head, too.

Item Number	Description	Pkg
G98035	Handheld Shower	4/cs

Replacement Tips

Heavy-duty rubber tip is reinforced with metal insert for longer wear. Fits walkers, commodes and transfer benches.

Item Number	Desc.	Size: Inch (cm)	Pkg
Replacement Tips			
MDS89755T	Gray	1 (2.5)	4 ea/st
Utility Tips			
G01182	Gray	3/4 (2)	8 set/cs
G01202	Gray	1 (2.5)	8 set/cs
G01212	Gray	1-1/8 (2.9)	8 set/cs

Item Number	Description	Color	Wt. Cap. Lbs (kg)	Seat Height Inch (cm)	Seat Depth Inch (cm)	Seat Width Inch (cm)	Pkg
Premium Padded Bath Bench							
G98027	with Back	Off-White	300 (136)	18-22 (46-59)	16.75 (42.5)	17 (43)	2 ea/cs
Standard Composite Bath Benches							
MDS89750R	without Back	White	300 (136)	17-21.5 (43-54.6)	16.75 (42.5)	22 (55.9)	2 ea/cs
MDS89755R	with Back	White	300 (136)	17-21.5 (43-54.6)	16.75 (42.5)	22 (55.9)	2 ea/cs
Basic Bath Bench							
G30403	without Back	White	250 (113)	13-20 (30-51)	12 (30.5)	19 (48)	4 ea/cs
G30402	with Back	White	250 (113)	13-20 (30-51)	12 (30.5)	19 (48)	3 ea/cs
G30401-4	without Back	Gray	250 (113)	13-20 (30-51)	12 (30.5)	19 (48)	4 ea/cs
G30400-4	with Back	Gray	250 (113)	13-20 (30-51)	12 (30.5)	19 (48)	4 ea/cs
MDS89740RW	without Back	Gray	250 (113)	14.5-19.5 (37-49.5)	12 (30.5)	20 (51)	2 ea/cs
MDS89745RA*	with Arms and Back	White	350 (159)	14.5-19.5 (37-49.5)	12 (30.5)	20 (51)	2 ea/cs
MDS89740RA	with Arms No Back	White	350 (159)	14.5-19.5 (37-49.5)	12 (30.5)	20 (51)	2 ea/cs
MDS89745R	with Back	Gray	250 (113)	14.5-19.5 (37-49.5)	12 (30.5)	20 (51)	2 ea/cs
MDS89740A	without Back	White	250 (113)	14.5-19.5 (37-49.5)	12 (30.5)	20 (51)	2 ea/cs
MDS89745A	with Back	White	250 (113)	14.5-19.5 (37-49.5)	12 (30.5)	20 (51)	2 ea/cs
Bariatric Bath Benches							
MDS89740AXW	without Back	White	550 (249)	14-17 (35.5-43)	12 (30.5)	20 (51)	Each
MDS89745AXW	with Back	White	550 (249)	16.5-20.5 (42-52)	21.5 (54.6)	18 (46)	Each

* Coming Soon

Retail Packaging Ready for Display

Padded Transfer Benches

- Allows patients to enter and exit tubs safely
- Padded seat ideal for long-term use
- Suction cup feet grip the tub floor for added safety
- Side arm provides extra stability and leverage



G98013A



G98338F



MDS86945P

Premium Padded Transfer Bench with Commode Opening

Commode opening affords the user greater comfort and better access for cleaning. Thick, high density foam cushions and backrest for comfortable seating that won't bottom out.

Standard Folding Padded Transfer Bench

Folding legs for easier transport and storage.

Basic Padded Transfer Bench

Comfortable, cushioned seat and back are covered in stain-resistant, light gray vinyl.

Item Number	Description	Color	Weight Cap Lbs (kg)	Seat Height Inch (cm)	Seat Width Inch (cm)	Seat Depth Inch (cm)	Pkg
G98013A	Padded with Commode Opening	Off-White	300 (136)	17-22 (43.2-55.9)	26.5 (67.3)	16.5 (41.9)	Each
G98338F	Folding Padded	Off-White	300 (136)	17-22 (43.2-55.9)	26.5 (67.3)	16.5 (41.9)	2/cs
MDS86945P	Basic Padded	Gray	300 (136)	20.5-25.5 (52.1-64.8)	16 (40.6)	17.25 (43.8)	2/cs

Transfer Bench Replacement Parts

Item Number	Description	Pkg
MDS86952LSC	Suction Cups with Leg Extensions	Pair
G01212	Utility Tips	8 st/cs

MDS86952LSC



G01212

Non-Padded Transfer Benches

GUARDIAN®

- Allows patients to enter and exit tubs safely
- Suction cup feet grip the tub floor for added safety
- Side arm provides extra stability and leverage
- Height-adjustable legs accommodate most tub wall heights



G98309G

Premium Transfer Bench

Built-in soap dish and hand-held shower slot maximize bathing convenience. Backrest adjusts without tools for right or left facing tubs. Comes in an attractive retail box.



MDS86960R

Standard Transfer Bench

No tools are required to assemble this durable plastic model that is rustproof. Included with this transfer bench is an optional removable arm that helps with patient balance and a built-in handheld shower holder and soap dish.



G98308A

Basic Transfer Bench

Textured, molded backrest and seating surface offer a non-slip surface with holes for water drainage. Backrest attaches for right or left facing tubs easily without tools.

Premium Bariatric Transfer Bench

Total of 6 legs and cross brace supporting legs inside tub provide maximum support for larger patients. Direction of the seat back can be reversed easily without tools.





MDS86952XW



MDS86952

Basic Transfer Bench

Blow molded transfer bench provides ideal support and stability at a more economical price. Direction of the seat back can be reversed easily without tools.

Item Number	Description	Color	Weight Cap Lbs (kg)	Seat Height Inch (cm)	Seat Width Inch (cm)	Seat Depth Inch (cm)	Pkg
Non-Padded Transfer Bench							
G98309G	Premium Non-Padded	White	300 (136)	17-21.5 (43.2-54.6)	31.5 (80)	21 (53.3)	3/cs 
MDS86960R	Standard Composite Plastic	White	350 (159)	17-21 (43.2-53)	17.5 (44.5)	29.5 (75)	2/cs 
G98308A	Non-Padded	Gray	300 (136)	15.5-20.5 (39.4-52.1)	26.5 (67.3)	16.25 (41.3)	2/cs
MDS86952	Basic	White	300 (136)	16.5-21.5 (41.9-54.6)	17 (43.2)	16.5 (41.9)	2/cs
Non-Padded Bariatric Transfer Bench							
MDS86952XW	Premium Bariatric	White	550 (250)	16-20 (40.6-50.8)	17 (43.2)	16 (40.6)	Each

 **Retail Packaging** Ready for Display

Elevated Toilet Seats

- Makes using the toilet easier for those with trouble bending
- Available in a variety of locking mechanisms and mounting methods



MDS80326

Extra Wide with Legs 600 lb (272 kg) Capacity

Steel legs go to the floor and soft foam armrests provide maximum support in getting on and off the toilet.



G30270A

Locking with Padded Arms

Padded arm rests assist in getting on and off toilet. Locking device secures seat to toilet and requires no tools.



MDS80320S

Locking with Molded Arms

Arms assist in getting on and off the seat. Locking device secures the seat to the toilet.



MDS80315

Universal Locking

Height-adjustable brackets allow seat to fit over any toilet. Brackets are vinyl covered to prevent chipping or scratching of porcelain. Seat is height-adjustable with plastic splash guard to make cleaning easy.



MDS80324S

Hinged Bolt-On








Seat can be lifted and lowered with hinged design. Existing seat and lid make this model less noticeable.



G30250

Basic Raised Seat

Simple design made of polypropylene resin is easy to clean. Non-slip pads and lip inside toilet offer additional stability.

Item Number	Description	Wt Capacity Lbs (kg)	Seat Elevation Inch (cm)	Pkg
With Arms				
MDS80326	Extra Wide with Legs	600 (272)	4 (10)	Each 
G30270A	Locking with Padded Arms	350 (159)	5 (12)	3/cs 
MDS80320S	Locking with Molded Arms	300 (136)	5 (12)	Each 
Without Arms				
MDS80315	Universal Locking	250 (113)	4.75-6.75 (12-17)	Each
MDS80324S	Standard Hinged Bolt-On	300 (136)	4 (10)	Each 
MDS80324E	Elongated Hinged Bolt-On	300 (136)	4 (10)	Each 
G30260*	Locking Raised	350 (159)	5 (12)	3/cs 
G30250	Basic Raised Seat	300 (113)	5 (12)	3/cs 

* Not Pictured

 **Retail Packaging** Ready for Display

Toilet Safety Frame



G30300

Item Number	Wt Capacity Lbs (kg)	Height x Width Inch (cm)	Pkg
G30300	250 (113)	26-31 (66-78.7) x 18-24 (46-61)	2/cs
G30300-1	250 (113)	26-31 (66-78.7) x 18-24 (46-61)	2/cs

Wall Grab Bars

GUARDIAN®

- Locate grab bars at any angle in the bathroom to help prevent falls
- Use grab bars in the kitchen, along stairways, or in the laundry room for additional fall prevention
- Combine grab bars with an elevated toilet seat or transfer bench for additional stability



Premium Knurled Chrome

Adjustable oval flange rotates to allow for secure mounting. Knurled finish provides steady grip even when wet.



Standard Knurled Chrome

Knurled-texture steel bars are easy to grip even in wet bathing areas.



Basic White Finish

White painted steel with white screws blend in to match most bathroom decors.

Item Number	Length Inch (cm)	Diameter Inch (cm)	Pkg
G81020KE	12 (30.5)	1 (2.5)	3/cs
G81030KE	16 (40.6)	1 (2.5)	3/cs
G81040KE	18 (45.7)	1 (2.5)	3/cs
G81050KE	24 (61.0)	1 (2.5)	3/cs
G81060KE	32 (81.3)	1 (2.5)	3/cs
G81070KE	36 (91.4)	1 (2.5)	3/cs

Item Number	Length Inch (cm)	Diameter Inch (cm)	Pkg
MDS86012CHRG	12 (30.5)	1 (2.5)	3/cs
MDS86016CHRG	16 (40.6)	1 (2.5)	3/cs
MDS86018CHRG	18 (45.7)	1 (2.5)	3/cs
MDS86024CHRG	24 (61.0)	1 (2.5)	3/cs
MDS86032CHRG	32 (81.3)	1 (2.5)	3/cs

Item Number	Length Inch (cm)	Diameter Inch (cm)	Pkg
G87520	12 (30.5)	1 (2.5)	3/cs
G87530	16 (40.6)	1 (2.5)	3/cs
G87540	18 (45.7)	1 (2.5)	3/cs
G87550	24 (61.0)	1 (2.5)	3/cs
G87560	32 (81.3)	1 (2.5)	3/cs
G87570	36 (91.4)	1 (2.5)	3/cs

Retail Packaging Ready for Display

Bath Tub Grab Bars and Safety Pole

- Give patients confidence and support with an extra grip to help them transfer into the tub



G98006

Deluxe Tub Assist Bar

Step through clamp design provides unobstructed entry/exit path allowing better hand positioning. Tool-free installation has stainless steel components with rubber-lined pads that won't mar tub surface. Adjust to fit tubs 2.75" - 6.25" (7-16 cm).



G98080

Bathtub Safety Rail

Clamps on two sides with foam lined clamp that won't mar tub surface. Mounting brackets adjust from 3.5" - 6" (9-15.24 cm) to fit most standard size tubs.



Floor to ceiling pole that safely assists standing, sitting, climbing stairs, or transferring.

G98120P

Item Number	Description	Pkg
G98006	Deluxe Tub Assist Bar	3/cs
G98080	Bathtub Safety Rail	3/cs
G98120B	Safe-T-Pole, Base	Each
G98120P	Safe-T-Pole, Pole	Each

Retail Packaging Ready for Display

Non-Folding Commodes

- Commodes can be used bedside for those who may not be mobile enough to walk to the bathroom
- Can also be used in the bathroom as an elevated toilet seat or as toilet safety rails when combined with most standard toilets
- All models include pail, lid and splash guard
- Clip-on seats can be removed easily for cleaning
- Gray powder-coated steel commodes are rust-resistant and durable
- Aluminum commodes are light weight and rustproof



G30211-4

Premium Chrome Plated Steel



MDS89664A

Standard Aluminum 3 in 1



G30213-4

Basic Painted Steel



Disposable Commode Liner

MDS89664LINER

Never clean a commode bucket again – for less than the price of a Sunday newspaper. Flip to page 64 for more information.



MDS89664

Basic Steel 3 in 1

Item Number	Description	Weight Capacity Lbs (kg)	Seat Height Range Inch (cm)	Seat Depth Inch (cm)	Width Between Arms Inch (cm)	Overall Width Inch (cm)	Pkg
G30211-4	Chrome Plated Steel	350 (159)	16.5-22.5 (41.9-57.2)	16 (40.6)	19 (48.3)	21.25 (54)	4/cs
MDS89664A	Aluminum	300 (136)	16-20 (41-51)	15 (38.1)	18.5 (47)	25.5 (65)	4/cs
G30213-4	Painted Steel	350 (159)	16.5-22.5 (41.9-57.2)	16 (40.6)	19 (48.3)	21.25 (54)	4/cs
MDS89664	Basic Steel	350 (159)	15-21 (38.1-53.3)	15 (38.1)	19 (48.3)	24 (61)	4/cs

Folding Commodes

GUARDIAN®

- Folding frame design assembles in seconds without tools
- Commode folds down to reduce storage space by 35% and makes delivery easier
- All models include pail, lid and splash guard
- Clip-on seats can be removed for easy cleaning



G30210-4F



MDS89664EL



G30213-4F

Premium Aluminum

Lightweight, aluminum design with folding capabilities makes this commode easy to bring in and out of storage.

Standard Elongated Seat

Features deeper seat depth to allow more room for cleaning.

Basic Steel

Standard folding commode helps save space at a more economic cost.

Item Number	Description	Weight Capacity Lbs (kg)	Seat Height Range Inch (cm)	Seat Depth Inch (cm)	Width Between Arms Inch (cm)	Overall Width Inch (cm)	Pkg
G30210-4F	Aluminum	250 (113)	17.5-21.5 (44.5-54.6)	16 (40.6)	19 (48.3)	21.25 (54)	4/cs
MDS89664EL	Elongated Seat	400 (182)	17-23 (43.2-58.4)	16.5 (41.9)	19 (48.3)	22 (55.9)	4/cs
G30213-4F	Steel	350 (159)	16.5-22.5 (41.9-57.2)	16 (40.6)	19 (48.3)	21.25 (54)	4/cs

Accessories

Item Number	Description	Pkg
MDS80306B	Replacement Commode Bucket	6/cs
MDS89664T	Commode Tips	4/cs
MDS89664SL	Seat and Lid Commode MDS89664	4/cs
MDS89664SLXW	Seat and Lid Commode MDS89664XW	Each
MDS89664SGC	Guard, Splash, Commode	4/cs



MDS80306B

Drop-Arm Commodes

- Arms swing out of the way for easier lateral transfers on and off the commode
- Transfers are easier for patients using a wheelchair
- Place near the bed with near arm down for patients transferring during the night
- Comes standard with pail, lids and splash guard



G98204

Premium Padded

Thick resilient foam padded seat and back for increased comfort is ideal for long-term use. Hard, chrome-plated steel is rust resistant and offers great durability.



G98202

Standard Non-Padded

Hard, chrome-plated steel frame is rust resistant and offers greater durability.



MDS89668

Basic Economy

Chrome-plated steel frame with easy, single-handed release mechanism for quick transfers. Comes standard with 12-quart pail with cover included.

Item Number	Description	Weight Cap Lbs (kg)	Seat Height Range Inch (cm)	Seat Depth Inch (cm)	Width Between Arms Inch (cm)	Overall Width Inch (cm)	Pkg
G98204	Premium Padded	350 (159)	20-25 (50.8-63.5)	18-20 (45.7-50.8)	18.5 (47)	22 (55.9)	Each
G98202	Standard Non-Padded	300 (136)	18.25-23.25 (46.4-59.1)	18-20 (45.7-50.8)	18.5 (47)	22 (55.9)	Each
MDS89668	Basic Economy	250 (113)	17-21 (43.2-53.3)	15 (38.1)	17.5 (44.5)	25.5 (64.8)	Each

Bariatric Commodes

GUARDIAN®

- Extra-wide seating area and higher weight capacities are ideal for larger patients
- Reinforced frames for greater weight-bearing capabilities



MDS89668XW

Bariatric Drop-Arm Commode 850 lb (386 kg) Capacity

Easy-to-release drop-arm mechanism allows larger patients to laterally transfer safely and with ease. Powder coated steel is rust resistant.



MDS89664XW

Bariatric Rigid Arm Commode 650 lb (295 kg) Capacity

Smaller frame, with high weight capacity takes up less space and fits over most standard toilets.



G30216B

Platform Bariatric Bedside Commode 500 lb (227 kg) Capacity

Platform seat and backrest provide additional comfort and support.



G30214-2

Heavy-Duty Commode 400 lb (181 kg) Capacity


Large base depth provides more stability.

Item Number	Description	Weight Cap Lbs (kg)	Seat Height Range Inch (cm)	Seat Depth Inch (cm)	Width Between Arms Inch (cm)	Overall Width Inch (cm)	Pkg
MDS89668XW	Bariatric Drop Arm	850 (386)	19.25 (48.9)	19 (48.3)	28 (71.1)	29 (73.7)	Each
MDS89664XW	Bariatric Rigid Arm	650 (295)	18.25-23.5 (46.4-59.7)	15 (38.1)	24 (61)	28 (71.1)	Each
G30216B	Platform	500 (227)	16.5-23.5 (41.91-59.7)	18.5 (47)	23.5 (59.7)	27 (68.6)	2/cs
G30214-2	Heavy Duty	400 (181)	17.5-21.5 (44.5-54.6)	16 (40.6)	21-1/8 (53.8)	23.75 (60.3)	2/cs

Never clean a commode bucket again—for less than the price of a Sunday newspaper.

Medline's disposable commode liners make commode cleanup easy and hygienic. Each liner includes an absorbent pad made from a polymer gelling agent that solidifies and neutralizes waste for convenient disposal. Because they are solidified, liquids will not move around in the commode, making transport easier. All you need to do is toss the liner in the garbage after use.

Disposable commode liners also eliminate contact between staff and human waste, to help reduce contamination. The average person will use two liners each day, at a suggest retail price of \$1.50 each.



Item Number	Description	Pkg
Disposable Commode Liner		
MDS89664LINER	Commode Liner	12 liners/box, 6 boxes/cs


How to use Disposable Commode Liners



1. Place the commode liner in the commode. Make sure that the liner is folded over the edge of the bucket.



2. As liquid hits the absorbent core, the polymer inside the core immediately begins converting the liquid into a solid. The core neutralizes odor as the liquid solidifies.



3. Simply toss the liner into the garbage after use.



"My goal is to bring the best possible service to my customers — keeping them profitable and helping them grow."

Kate M., Medline Sales Representative



Medline Industries, Inc.

One Medline Place
Mundelein, IL 60060

Medline United States **1-800-MEDLINE (633-5463)**

www.medline.com | info@medline.com

Medline Canada **1-800-396-6996**

www.medline.ca | canada@medline.com

Medline México **01-800-831-0898**

www.medlinemexico.com | mexico@medline.com

Some products may not be available for sale in Mexico or Canada. Please contact your Account Representative for additional information. We reserve the right to correct any errors that may occur within this catalog. ©2013 Medline Industries, Inc. Medline, Easy Care, Biomask and Simplicity are trademarks and Red Dot, Excel Translator, Excel Hybrid 2, CURAD, Guardian, Excel, Accu-Therm, Therahoney, Fitright, Remedy, and Quick Alert are registered trademarks of Medline Industries, Inc. ArmaSport is a registered trademark of Armacell Enterprise GmbH Corp. Rifton is a registered trademark of Community Products, LLC, Linak is a registered trademark of Linak A/S LLC, Koolnit is a registered trademark of M.W. Sales and Service Corp. Microban is a registered trademark of Microban Products Company Corp. Isoflo Memory Control and Quadtro Select are registered trademarks of Roho, Inc. InTENSity is a registered trademark of Roscoe Medical, Inc. UMP is a registered trademark of Stanley Security Solutions, Inc. Quadra3D is a registered trademark of The Comfort Companies, Inc. Tinnerman is a registered trademark of Trans Technology Engineered Components, Nylex is a registered trademark of Xamax Industries, Inc. Leg Loop is a trademark of Maddak. Coldspot is a trademark and The Cuff, Point Relief, Hotspot, Digi-Extend, Theraputty, Digi-Flex, Digi-Squeeze are registered trademarks of Elliott Goldberg.

MKT1331882 / LIT124R / 10M / QG / 10