

Point of Care

Testing Solutions

A^{Pro}
ADVANTAGE[®]

by **NDC**

Pregnancy

Rapid, one-step test for the qualitative detection of human chorionic gonadotropin (hCG).

Infectious Diseases

A rapid chromatographic immunoassay test for the qualitative detection of Strep A antigen from throat swab specimens. A rapid chromatographic immunoassay technology identifies heterophile antibodies in whole blood, serum or plasma to aid in the diagnosis of Infectious Mononucleosis.

Cancer Screening

Convenient and qualitative methods for detection of asymptomatic occult bleeding. Available in guaiac slide.

Gastrointestinal Disorders

A rapid test for the qualitative detection of IgG antibodies to *Helicobacter pylori* in whole blood, serum and plasma.



Pregnancy Test Devices

Pro
ADVANTAGE[®]

by **NDC**

Rapid, one step tests for the qualitative detection of human chorionic gonadotropin (hCG) in serum and or urine.

Three variations of Pro Advantage hCG tests are available:
Urine Cassette Devices, Urine/Serum Pregnancy Devices and Urine Test Strips.

- >99% Sensitivity; >99% Specificity
- 25 mIU/mL minimum hCG detection level
- Accurate results in 3 minutes for urine and 5 minutes for serum
- CLIA waived for urine and FDA 510(k) cleared to market.



Item #	Description
--------	-------------

Pregnancy Test Devices

P080063	Urine hCG Pregnancy Cassette Device
----------------	--

Includes: 25 individually packaged devices and droppers. 25 tests per box

P080021	Urine/Serum hCG Pregnancy Cassette Device
----------------	--

Includes: 25 individually packaged devices and droppers. 25 tests per box

P080023	Urine hCG Pregnancy Test Strips
----------------	--

Includes: 20 Test Strips per canister, 5 canisters per box



**For Urine Samples,
Serum:
Moderately Complex**
CPT Code:
81025 Urine, 84703 Serum

Strep A Test



A rapid test for the qualitative detection of Strep A antigen in throat swab specimens.

Lateral flow immunoassay technology identifies antibodies specific for whole cell Lancefield Group A Streptococcus to selectively detect Strep A carbohydrate antigen in a throat swab specimen. Aids in the early diagnosis of Group A Streptococcal infection.

- 97% Sensitivity; 95% Specificity
- Accurate results in 5 minutes
- Quickly diagnose and treat infections to help avoid serious complications
- CLIA waived and FDA 510(k) cleared to market
- Room temperature storage



Item #	Description
--------	-------------

P080033	Strep A Test-Throat Swab Kit <i>Includes: 25 test strips, 25 swabs, 25 extraction tubes, 1-10mL reagent A, 1-10mL reagent B, 1-1mL positive control, 1-1mL negative control, package insert, procedure card. 25 tests per box</i>
----------------	---



CPT Code: 87880

Infectious Mononucleosis Test Device



Rapid chromatographic immunoassay technology identifies heterophile antibodies in whole blood, serum or plasma to aid in the diagnosis of Infectious Mononucleosis.

- Qualitative membrane strip based immunoassay for detection of IM heterophile antibodies
- All inclusive 20 test kit
- Simple to perform using whole blood (from venipuncture or fingerstick), serum or plasma
- Reliable results in 5 minutes
- CLIA waived for whole blood and FDA 510(k) cleared to market



Item #	Description
--------	-------------

P080016	Infectious Mononucleosis Test Device <i>Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box</i>
----------------	--



For Whole Blood
Serum/Plasma:
Moderately Complex
CPT Code: 86308

Fecal Occult Blood Test Devices

A^{Pro} ADVANTAGE[®]

by **NDC**

Rapid, convenient, and qualitative methods for asymptomatic occult bleeding. Available in guaiac slide – with both color enhanced sensitivity and readability.

- On board positive and negative monitors
- Includes enhanced color sensitivity and readability
- CLIA waived and FDA 510(k) cleared to market



Item #	Description
--------	-------------

Fecal Occult Blood Test Devices

P080008 Fecal Occult Blood Patient Pack – Guaiac Slide

Includes: 100 Triple Slide Tests, 6X10mL Developers, Applicators, Foil-Lined U.S. Postal Service approved envelopes. 2x50 per box

P080018 Fecal Occult Blood Lab Pack – Guaiac Slides

Includes: 100 Single Slides Tests, 2x10mL Developers, Applicators. 100 per box

P080028 Fecal Occult Developer Solution – for Guaiac Slides

Includes: 8x10mL Developers. 8 per box



CPT Code: 82270

H. pylori Test Device

A^{Pro} ADVANTAGE[®]

by **NDC**

A one step, qualitative membrane strip based immunoassay for the detection of H. pylori IgG antibodies in whole blood, serum or plasma.

H. pylori is one of the most common bacterial infections in humans causing active and chronic gastritis, non-ulcer dyspepsia and gastric/duodenal ulcer formation. The Pro Advantage H. pylori Test Device – Whole Blood/Serum/Plasma offers a simple and quick test format which assists the physician in making a treatment decision during the same office visit.

- 88% sensitivity; 89% specificity
- Accurate results in 10 minutes
- Reduces the risk of developing more serious illnesses
- Lessens the likelihood of complications
- CLIA waived for whole blood and FDA 510(k) cleared to market
- Fast and easy – no biopsy or histology expertise required



Item #	Description
--------	-------------

P080095 H. pylori Test Device – Whole Blood/Serum/Plasma

Includes: 20 individually packaged devices, 20 disposable sample droppers, 20 disposable heparinized capillary tubes, 1 dispensing bulb, 1 sample buffer, 1 procedure card, 1 package insert. 20 tests per box



**For Whole Blood,
Serum/Plasma:
Moderately Complex**

CPT Code: 86677

Distributed By: