

Steel Frame Privacy Screen

3130/3420 Standard White

3139/3429 Sure-Chek®



ISO 9001-2008 Certified

Create a NON-CLINICAL ENVIRONMENT for your patients with Winco's Sure-Chek® colored vinyl. Style, Versatility, & Durability!

Winco's steel frame privacy screens stand up to years of hard use. The wide, sturdy and durable steel frame helps to eliminate tipping or leaning. Both screen models come with twin wheel casters that can be easily converted to crutch tip bases.



KEY FEATURES

- Wide steel frame
- Twin casters
- Replaceable screen panels

1.800.237.3377 • 5516 SW 1st Lane • Ocala, Florida 34474

www.wincomfg.com



PRODUCT SPECIFICATIONS

STEEL FRAME PRIVACY SCREEN

3130/3420 STANDARD WHITE

3139/3429 SURE-CHEK®


- 1 Choose from our bacteriostatic vinyl panels or our OPTIONAL Sure-Chek® antimicrobial material.
- 2 Sturdy hinged panel ensures years of problem-free use.
- 3 Wide base for added stability
- 4 Comes with smooth rolling twin wheel casters or use with optional crutch-tips.
- 5 Smooth powder-coated steel frame for durability and ease of cleaning.



SPECIFICATIONS

	3130/3139	3420/3429
Base Size	20" x 27" (50.8 x 68.9cm)	20 x 42" (50.8 x 106.7cm)
Overall Height	70" (177.8cm)	69" (175.3cm)
Overall Width	79" (200.7cm)	42" (106.7cm)
Single Panel Length	49" (124.46cm)	49.5" (125.73cm)
Single Panel Width	27" (68.6cm)	42" (106.7cm)
Weight	25 lbs. (11.3kg)	15" (6.8kg)

COLORS



Optional Optional Optional Optional Optional

Printed colors shown may vary slightly from actual color of product.
 Sure-Chek® vinyls are antimicrobial, antistatic, odor resistant, stain resistant, & fire retardant.